

Biomarker discovery: LC-MS Proteomics


Nanocourse: Data Science using R
September 6th, 2024

Jeon Lee

Agenda

1. Overview of biomarker discovery steps
2. Batch correction/data harmonization
3. Harmonization of proteomics data with missing values
4. Introduction to LC-MS proteomics
 - MS for metabolomics/proteomics
 - LC; LC-MS/MS
 - Peak annotation
 - Typical proteomics data & analysis steps
5. Demo: Proteomics analysis
6. Hands-on practice

Overview of biomarker discovery steps




<Adopted from https://app.biorender.com/profile/auriol_purdie/templates/6635bb7572a44e4ad29ce313>

For RNAseq, have all the data

Batch correction/data harmonization (1)

To generalize findings with high confidence, need more samples.




Batch effect- Comparing two samples is not enough due to batch differences

Yu, Y., Zhang, N., Mai, Y. *et al.* Correcting batch effects in large-scale multiomics studies using a reference-material-based ratio method. *Genome Biol* **24**, 201 (2023).

Batch correction/data harmonization (2)


a Transcriptomics



Has slightly better ability to remove batch effect (could be bc has more features/data points than proteomics (10000 vs 1000))

Use similar technique, but how batch effect influences each level depends


b Proteomics




Yu, Y., Zhang, N., Mai, Y. *et al.* Correcting batch effects in large-scale multiomics studies using a reference-material-based ratio method. *Genome Biol* **24**, 201 (2023).

Batch correction/data harmonization (3)

C Metabolomics



d





Yu, Y., Zhang, N., Mai, Y. *et al.* Correcting batch effects in large-scale multiomics studies using a reference-material-based ratio method. *Genome Biol* **24**, 201 (2023).

Usually, lots of missing data points

Harmonization of proteomics data with missing values (1)

For single cell and LC-MS data:
Have missing values

Can't batch correct imputed data bc need complete data to batch correct. So
imput first




Always have missing values


Voß, H., Schlumbohm, S., Barwikowski, P. *et al.* HarmonizR enables data harmonization across independent proteomic datasets with appropriate handling of missing values. *Nat Commun* **13**, 3523 (2022)

Harmonization of proteomics data with missing values (2)

Volcano plot visualization of individual experimental setups



Volcano plot comparison of different strategies




Voß, H., Schlumbohm, S., Barwikowski, P. *et al.* HarmonizR enables data harmonization across independent proteomic datasets with appropriate handling of missing values. *Nat Commun* **13**, 3523 (2022)

Introduction to LC-MS proteomics


MS for metabolomics/proteomics

- MS is an analytical technique that ionizes chemical species and sorts the ions based on their mass-to-charge ratio (m/z).
- Mass spectrometers are comprised of an ionization source and a mass detector, for example,
 - MALDI-TOF: matrix assisted laser desorption ionization, time-of-flight detection
 - ESI-trap: electrospray ionization, ion trap detection




Liquid chromatography

- LC is the separation technique of choice for larger and non-volatile molecules such as proteins and complex peptides
- LC is also an ideal method for separating isomers, which have the same mass and will otherwise not be differentiated by a mass spectrometer
- LC-MS offers broad sample coverage because different column chemistries, such as reversed phase liquid chromatography, can be used




Ex. LC-MS chromatogram for 32 targeted analytes
Swortwood, M. et al. (2012), Anal. Bioanal. Chem.



Tandem mass spectrometry (LC-MS/MS)



- Combination of LC and two mass analyzers in mass spectrometry (MS/MS)
- Once samples are ionized to generate a mixture of ions, precursor ions of a specific mass-to-charge ratio (m/z) are selected (MS1) and then fragmented (MS2) to generate a product ions for detection.
- The fragments then reveal aspects of the chemical structure of the precursor ion.



To get better res, ionizes twice

Tandem_mass_spectrometry, Wikipedia

A brief history of mass spectrometry,
healthcare-in-europe.com




F.A. Mellon, Encyclopedia of Food Sciences and
Nutrition (2nd Edition, 2003)

UT Southwestern
Medical Center
Lyda Hill Department of Bioinformatics

CTTT lab

Peak annotation

- Peak grouping (a-c) aims at grouping peaks that belong to each metabolite/protein.
- In feature annotation (e, f), expected theoretical distances between known ion adduct masses are compared with experimental distances found among peaks (e).
- After peak annotation, putative identification can be achieved by accurate mass search (g) or by comparison with MS/MS data (h).



Xavier, D. et al (2018) Analytical chemistry


Typical proteomics data

- Each row includes an accession id, of which protein has been detected, other meta data, and its abundance measured across the samples.
 - Coverage[%] = no. amino acids in all found peptides / total no. amino acids in the entire protein sequence
 - No.Peptides: No. distinct peptide sequences in the protein group
 - No.Unique_Peptides: No. peptide sequences unique to a protein group
 - No.PSMs: Total no. identified peptide sequences for the protein, including those redundantly identified
 - MW[kDa]: Molecular weight without considering post-translational modifications
- Some proteins are detected in some samples but not in the other samples

| | A | C | D | E | F | G | H | I | J | K | L |
|----|-----------|-----------|--|-------------|-------------|---------|--------------------|---------|--------|------------------|------------------|
| 1 | Protein_F | Accession | Description | Coverage[%] | No.Peptides | No.PSMs | No.Unique_Peptides | MW[kDa] | Gene | Abundance_934187 | Abundance_934188 |
| 2 | High | P63261 | Actin, cytoplasmic 2 OS=Homo sapiens OX=9606 GN=ACTG1 PE=1 SV=1 | 96 | 43 | 6502 | 2 | 41.8 | ACTG1 | 30142523.02 | 449237770.5 |
| 3 | High | P60709 | Actin, cytoplasmic 1 OS=Homo sapiens OX=9606 GN=ACTB PE=1 SV=1 | 96 | 43 | 6488 | 2 | 41.7 | ACTB | 11088943097 | 1.28E+11 |
| 4 | High | O43707 | Alpha-actinin-4 OS=Homo sapiens OX=9606 GN=ACTN4 PE=1 SV=2 | 95 | 113 | 5193 | 6 | 104.8 | ACTN4 | 2053379189 | 54971240930 |
| 5 | High | Q13813 | Spectrin alpha chain, non-erythrocytic 1 OS=Homo sapiens OX=9606 GN=SP | 89 | 319 | 5110 | 12 | 284.4 | SPTAN1 | 3401878545 | 26753571901 |
| 6 | High | A0A0D9SF | Spectrin alpha chain, non-erythrocytic 1 OS=Homo sapiens OX=9606 GN=SP | 88 | 309 | 4996 | 2 | 282.7 | SPTAN1 | 6749205.641 | 42286992 |
| 7 | High | H7C144 | Alpha-actinin-4 OS=Homo sapiens OX=9606 GN=ACTN4 PE=1 SV=2 | 92 | 111 | 4830 | 5 | 104.3 | | 1874029 | 52463477.09 |
| 8 | High | Q01082 | Spectrin beta chain, non-erythrocytic 1 OS=Homo sapiens OX=9606 GN=SPT | 84 | 246 | 4266 | 224 | 274.4 | SPTBN1 | 2552146613 | 20243849937 |
| 9 | High | P68032 | Actin, alpha cardiac muscle 1 OS=Homo sapiens OX=9606 GN=ACTC1 PE=1 SV | 64 | 32 | 3429 | 9 | 42 | ACTC1 | 906591877.7 | 9311096952 |
| 10 | High | P12814 | Alpha-actinin-1 OS=Homo sapiens OX=9606 GN=ACTN1 PE=1 SV=2 | 88 | 82 | 2375 | 3 | 103 | ACTN1 | 155772807.2 | 5159761723 |
| 11 | High | A0A7I2V4 | Alpha-actinin-1 OS=Homo sapiens OX=9606 GN=ACTN1 PE=1 SV=1 | 83 | 77 | 2313 | 0 | 103.5 | | | 1527293.125 |
| 12 | High | Q15149 | Plectin OS=Homo sapiens OX=9606 GN=PLEC PE=1 SV=3 | 71 | 375 | 2167 | 131 | 531.5 | PLEC | 1094893182 | 4176426088 |

Proteomics analysis steps

- Three steps in proteomics analysis
 - Data acquisition/raw data processing
 - Data pre-processing
 - (main/downstream) data analysis
- Typical data pre-processing steps include: 1) missing value imputation, 2) transformation, (3) scaling, and (4) normalization



Demo: Proteomics analysis

“Demo_Proteomics_Analysis.html”